



## CD45 RA/RO Panel (SCID)

Order Name: **CD45 RA/RO**  
Test Number: 3654613  
Revision Date: 02/19/2024

| TEST NAME               | METHODOLOGY    | LOINC CODE |
|-------------------------|----------------|------------|
| CD45 RA/RO Panel (SCID) | Flow cytometry |            |

### SPECIMEN REQUIREMENTS

|              |   |                    |                            |                         |
|--------------|---|--------------------|----------------------------|-------------------------|
| Specimen     | Specimen Volume (min)   | Specimen Type      | Specimen Container         | Transport Environment   |
| Preferred    | <b>4 mL (0.5 mL)</b>  | <b>Whole Blood</b> | <b>EDTA (Lavender Top)</b> | <b>Room Temperature</b> |
| Instructions | <p>Collect Whole blood <b>EDTA (Lavender Top)</b>.<br/>Keep whole blood, do not centrifuge or refrigerate.<br/><b>Deliver to the laboratory ASAP - Same Day as Collection!</b><br/><b>Collect Monday through Thursday only. Specimens must be received at Labcorp Oklahoma Inc. by Thursday afternoon or collect on Monday.</b></p> |                    |                            |                         |

### GENERAL INFORMATION

|                     |  |
|---------------------|--|
| Testing Schedule    | Tue-Sat  |
| Expected TAT        | 2 days   |
| Notes               | This panel is also known as: Lymphocyte Subset 6 - Total Lymphocyte Enumeration with CD45RA and CD45RO               |
| CPT Code(s)         | 86355, 86357, 86359, 86360, 86356x2  |
| Service Provided By |  <b>labcorp</b><br>Oklahoma, Inc. |