


Volatiles, Whole Blood

Order Name: **Volatiles WB**
Test Number: 4302200
Revision Date: 07/07/2025

TEST NAME	METHODOLOGY	LOINC CODE
Acetone Level	<u>Gas chromatography (GC)</u>	5568-1
Methanol Level	<u>Gas chromatography (GC)</u>	5693-7
Isopropanol Level	<u>Gas chromatography (GC)</u>	5669-7
Ethanol Level	<u>Gas chromatography (GC)</u>	56478-1

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	7mL (0.5mL)	Whole Blood	Sodium Fluoride (Gray)	Room Temperature
Alternate 1	7mL (0.5mL)	Serum	No Additive Clot (Red Top, No-Gel, Plastic)	Room Temperature
Instructions	<p>Specimen: 7 mL (0.5 mL) Whole blood or serum - Room Temperature</p> <p>Container: Gray-top (sodium fluoride) tube (preferred) or red-top tube. Submit original unopened tube.</p> <p>Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 14 days (Freeze/thaw cycles Stable x3)</p> <p>Special Instructions: This test number is not intended for workplace testing and does not comply with state regulatory workplace testing programs.</p> <p>Do not prepare venipuncture site with alcohol or remove stopper from tube.</p>			

GENERAL INFORMATION	
Expected TAT	2-3 days
CPT Code(s)	80320
Service Provided By	 labcorp Oklahoma, Inc.