


Pseudoephedrine Confirmation Urine

Order Name: **PSEUDO CNF**
Test Number: 4000138
Revision Date: 05/14/2020

TEST NAME	METHODOLOGY	LOINC CODE
Pseudoephedrine Confirmation	Liquid Chromatography/Mass Spectroscopy	4004-8
Pseudoephedrine Interpretation	Interpretive information.	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, Random	Sterile Container or Monovette	Refrigerated
Instructions	Keep refrigerated. Specimen Stability: Ambient 7 days. Refrigerated 30 days. Frozen 6 months Please provide accurate prescription list to lab.			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	3-5 days
Clinical Use	Used to evaluate medication compliance in pain management.
CPT Code(s)	G0480 (alternate:80324)
Service Provided By	 labcorp Oklahoma, Inc.