

AH50 (Complement Activity, Alternative Pathway)

Order Name: **AH50**
Test Number: 5002275
Revision Date: 04/14/2023

TEST NAME	METHODOLOGY	LOINC CODE
AH50 (Complement Activity, Alternative Pathway)	Hemolysis	13087-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p>Specimen Type: Red-top tube</p> <p>Specimen Storage: FREEZE (SHIP ON DRY ICE) Do Not Share this specimen with another test.</p> <p>Specimen Collection: Transfer specimen to plastic transport tube before freezing. Send to lab on dry ice.</p> <p>Special Instructions: Testing referred to National Jewish Health NJHCO#AM50</p> <p>Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION	
Expected TAT	27 - 36 days
Notes	Labcorp Test Code: 831782
CPT Code(s)	86161
Service Provided By	 labcorp Oklahoma, Inc.