

5-HIAA Qn,Random Urine

Order Name: **5-HIAA U Rand**
Test Number: 5194552
Revision Date: 12/12/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|-------------------------------|-------------|------------|
| 5-HIAA Qn,Random Urine | | |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|---------------------|---|---------------|--------------------|-----------------------|
| Instructions | Specimen Type: Plastic urine container, no preservative (NOTE: Boric acid may be added as a preservative for other tests without harm to 5-HIAA.) Specimen Storage: Room Temperature Specimen Collection: Screw the lid on securely. Label container with patient's name, date, and time. Mix well; send aliquot. Specimen Stability: Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days | | | |

GENERAL INFORMATION

| | |
|---------------------|---|
| Notes | Labcorp Test Code: 316205 |
| CPT Code(s) | 82570, 83497 |
| Service Provided By |  labcorp Oklahoma, Inc. |