

**6-Acetylmorphine, MS, Ur, MAT**

Order Name: **6 AM UR**  
Test Number: **5195572**  
Revision Date: **01/09/2024**

TEST NAME	METHODOLOGY	LOINC CODE
6-Acetylmorphine, MS, Ur, MAT	<a href="#">Mass Spectrometry (MS)</a>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.			

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80356
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.