

Drug Screen 17 w/Conf, UR

Order Name: **Drug 17 w/Conf, UR**
Test Number: 5195607
Revision Date: 01/09/2024

| TEST NAME | METHODOLOGY | LOINC CODE |
|---------------------------|------------------|------------|
| Drug Screen 17 w/Conf, UR | Immunoassay (IA) | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|--|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 30 mL (15 mL) | Urine, Random | Sterile Urine Container Preservative Free | Room Temperature |
| Instructions | Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks. | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Expected TAT | 5-8 days |
| CPT Code(s) | 80307 |
| Service Provided By |  labcorp Oklahoma, Inc. |